

The Commonwealth of Massachusetts
Executive Office of Health and Human Services
Department of Public Health
William A. Hinton State Laboratory Institute
305 South Street, Jamaica Plain, MA 02130

DEVAL L. PATRICK
GOVERNOR

TIMOTHY P. MURRAY
LIEUTENANT GOVERNOR

JUDYANN BIGBY, MD
SECRETARY

JOHN AUERBACH
COMMISSIONER

01/26/2012

Elle Rackemann
Assistant District Attorney, Suffolk County

Dear ADA Rackemann,

Enclosed is the information you requested in regards to Commonwealth vs. [REDACTED]
Included are copies of the following:

1. Curriculum Vitae for Annie Dookhan and Nicole Medina.
2. Drug Analysis Laboratory Receipt.
3. Control Cards with analytical results for samples # [REDACTED]
4. Analysis sheets with custodial chemist's hand notations and test results.
5. GC/Mass Spectral analytical data for samples # [REDACTED]

Annie Dookhan was the custodial chemist and performed the preliminary testing and net weight for this sample. Nicole Medina analyzed the GC/MS data for this sample.

If you have any questions about these materials, please call me at the number below.

Sincerely,

A handwritten signature in black ink, appearing to read "Annie Dookhan".

Annie Dookhan
Chemist
Drug Analysis Lab
Jamaica Plain, MA. 02130
(617) 983-6622

Curriculum Vitae

Annie Khan (Dookhan)

Education:

University of Massachusetts, Boston, Ma, Master of Science in Chemistry

University of Massachusetts, Boston, Ma, Bachelor of Science

Major: Biochemistry

Minor: Economics

Experience:

2003 – present

Chemist I, II, Massachusetts Department of Public Health, Drug Analysis Laboratory

*Responsible for the identification of illicit drugs to determine violations of harmful and narcotic drug laws.

*Writing of the Standard Operating Procedures (SOPs) and Protocols.

*Responsible for the Quality Control and Quality Assurance of analytical instruments, reagents and controls/standards.

*Maintenance and repairs of analytical instruments.

*Trained in the use of complex analytical instrumentation, microscopes and balances for the purpose of drug analysis.

*Qualified as an expert witness in Massachusetts Courts and U.S. District Court.

*Completed six-week training course conducted by senior staff within the Department of Public Health, Drug Analysis Laboratory.

*Appointed Assistant Analyst by Assistant Commissioner of Public Health.

*Notary Public.

2001 – 2003

QC Analyst I, II, UMMS-Massachusetts Biologic Laboratory, QC Material Control

*Completed proficiency training conducted by a member of the staff within the MBL Quality Control and Quality Assurance Department for the Food and Drug Administration.

*Method Development for creating new techniques for the QC Dept. and FDA.

*Writing, revising and reviewing the Standard Operating Procedures (SOPs).

*Trained and supervised new chemists and interns for the department.

*Routine QC testing of products for the FDA.

*Trained in the use of complex analytical instrumentation, and balances for the purpose of QC analysis for product and validation projects.

*Calibration, preventive maintenance, QC and QA of analytical instrumentation.

*Complete testing of chemicals for Vendor Validation Project for the FDA.

*Compendial testing and interpretation of the USP, ACS, FCC, AOAC, Merck Index, PDR, etc.

Additional Courses:

Department of Justice – Forensics Professionals

Introduction to Forensic Drug Chemistry

Chemical Spot Test and Mechanism for Illicit Drugs

Forensic Mass Spectrometry

Forensic Science 101: An Introduction

Fundamentals of Forensic Questioned Documents

Perspectives in Expert Testimony

Evidence Analysis

Ethics in Forensic Science

Transition to Leadership

Fundamentals of Forensic Toxicology

Additional Trainings:

Good Laboratory Practice/Good Manufacturing Practice (GLP/GMP) training with Massachusetts Biologic Laboratory.

Quality Control/Quality Assurance (QC/QA) training according to FDA Codes and Regulations.

Gas Chromatography (GC) and Gas Chromatography/Mass Spectrometry (GC/MS) trainings with Agilent Technologies and Restek.

High Power Liquid Chromatography (HPLC) and HPLC/MS/MS trainings with Waters Cooperation.

Fourier Transform Infrared Spectroscopy (FTIR) training with Spectros and Spectrum One.

Ultraviolet-Visible Spectroscopy (UV-Vis) training with Agilent Technologies and Beckman.

Total Organic Carbon Analyzer (TOC) training with MBL and Sievers.

Associations:

American Chemical Society (ACS)

Northeastern Association of Forensics Science (NEAFS)

Curriculum Vitae

Nicole E. Medina

Education

Bachelor of Science, Chemistry (1995)

Emmanuel College, Boston, MA

Research Project: (1994)

In-Vivo Mapping of Cocaine Metabolites in Biological Fluids

Conducted at Massachusetts State Police Crime Lab, Toxicology Department

Employment

Chemist II State Laboratory Institute (March 2008-Present)

Massachusetts Department of Public Health

Drug Analysis Laboratory

- Responsible for the identification of substances and trafficking substances to determine violation of the Massachusetts drug laws
- Responsible for the identification of pharmaceuticals to determine violation of the Massachusetts drug laws
- Operate analytical instrumentation, microscopes and balances for forensic drug analysis

Chemist I State Laboratory Institute (2004-March 2008)

Massachusetts Department of Public Health

Drug Analysis Laboratory

- Responsible for the identification of substance to determine violation of the Massachusetts drug laws
- Operate analytical instrumentation, microscopes and balances for forensic drug analysis
- Successfully completed an eight week training in the analysis of drugs conducted by senior staff of the Department of Public Health, Drug Analysis Laboratory
- Appointed assistant analyst for the Department of Public Health, Drug Analysis Laboratory in 2005.

Research Associate

(September 2001 – October 2004)

Wolfe Laboratories, Newton & Watertown, MA

Laboratory Technician

(January 1996 - April 1997)

Gillette, Boston, MA

The Commonwealth of Massachusetts

Executive Office of Health and Human Services

Department of Public Health

State Laboratory Institute

Boston Drug Laboratory

Tel (617) 983-6622

Fax (617) 983-6625

Amherst Drug Laboratory

Tel (413) 545-2601

Fax (413) 545-2608

Boston Hours

8:00 – 11:00

2:00 – 4:00

Amherst Hours

9:00 – 12:00

1:00 – 3:00

DRUG RECEIPT

City or Department: Boston

Police Reference No.: [REDACTED]

Name and Rank of Submitting Officer: PO Sybil White

11064

Defendant(s) Name (last, first, initial):

[REDACTED]

To be completed by Submitter
Description of Items Submitted

To be completed by Lab Personnel
Gross Weight Lab Number

5 pills

9.28 gm

5 pills

13.63 gm

4 pills

9.17 gm

39 pills

49.97 gm

AKL

RESUBMITTED TO BE TESTED

Received by: AS

Date: 7.13.11

The Commonwealth of Massachusetts
Executive Office of Health and Human Services
Department of Public Health
State Laboratory Institute

Boston Drug Laboratory
Tel (617) 983-6622
Fax (617) 983-6625

Boston Hours
8:00 – 11:00
2:00 – 4:00

Amherst Drug Laboratory
Tel (413) 545-2601
Fax (413) 545-2608

Amherst Hours
9:00 – 12:00
1:00 – 3:00

DRUG RECEIPT

City or Department: Boston Police Reference No.: [REDACTED]

Name and Rank of Submitting Officer: PO Stephen White 11064

Defendant(s) Name (last, first, initial):
[REDACTED]

To be completed by Submitter
Description of Items Submitted

To be completed by Lab Personnel
Gross Weight Lab Number

5 pills	9.28 gm	[REDACTED]
5 pills	13.63 gm	
4 pills	9.17 gm	
39 pills	49.97 gm	
ALL		
RESUBMITTED TO BE TESTED		

Received by: AS Date: 7.13.11

No. [REDACTED] Date Analyzed: 09/19/11
City: Boston D.C.U. Police Dept.
Officer: P.O. SYBIL WHITE
Def: [REDACTED]
Amount: 5.0 Subst: TABS
No. Cont: 1 Cont: pb
Date Rec'd: 07/13/2011 No. Analyzed:
Gross Wt.: 9.28 Net Weight:
Tests: 1 ASD
2DEM
Prelim: clonidine
Findings: Clonidine
E

No. [REDACTED] Date Analyzed: 09/19/11
City: Boston D.C.U. Police Dept.
Officer: P.O. SYBIL WHITE
Def: [REDACTED]
Amount: 5.0 Subst: TABS
No. Cont: 1 Cont: pb
Date Rec'd: 07/13/2011 No. Analyzed:
Gross Wt.: 13.63 Net Weight:
Tests: 1 ASD
2DEM
Prelim: gabapentin
Findings: Gabapentin
E

No. [REDACTED] Date Analyzed: 09/19/11
City: Boston D.C.U. Police Dept.
Officer: P.O. SYBIL WHITE
Def: [REDACTED]
Amount: 4.0 Subst: TABS
No. Cont: 1 Cont: pb
Date Rec'd: 07/13/2011 No. Analyzed:
Gross Wt.: 9.17 Net Weight:
Tests: 1 ASD
2DEM
Prelim: clonidine
Findings: Clonidine
E

No. [REDACTED] Date Analyzed: 09/19/11
City: Boston D.C.U. Police Dept.
Officer: P.O. SYBIL WHITE
Def: [REDACTED]
Amount: 39.0 Subst: TABS
No. Cont: 1 Cont: pb
Date Rec'd: 07/13/2011
Gross Wt.: 49.97
No. Analyzed:
Net Weight:
Tests: 1 ASD
2DEM

Prelim: gabapentin
Findings: Gabapentin
E

DRUG POWDER ANALYSIS FORM

SAMPLE # [REDACTED] AGENCY Boston ANALYST ASJ

No. of samples tested: Evidence Wt.

PHYSICAL DESCRIPTION:

5 white round tablet
in 1 pb

Gross Wt ():

Gross Wt ():

Pkg. Wt:

Net Wt: 0.4595



PRELIMINARY TESTS

Spot Tests

Cobalt
Thiocyanate ()

Marquis

Froehde's

Mecke's

Microcrystalline Tests

Gold
Chloride

TLTA ()

OTHER TESTS

Muonexor A+ (+)

PRELIMINARY TEST RESULTS

RESULTS clonidine

DATE 08-31-11

GC/MS CONFIRMATORY TEST

RESULTS clonidine

MS
OPERATOR DEM

DATE 9-19-11

Drug Information

DRUG NAME: CLONIDINE HYDROCHLORIDE

INGREDIENTS: CLONIDINE HYDROCHLORIDE -- 0.1 MG

Color: White

Shape: Circle

Imprint: MYLAN 152

Form: Oral Tablet

Available-Container-Size: Strip of 90 , Strip of 25 , Bottle of 100 , Bottle of 1000 , Strip of 1 , Strip of 30 , Strip of 100

AAPCC Code: 201081 - Clonidine

NDC: 51079-0299-62

51079-0299-19

00378-0152-01

00378-0152-10

51079-0299-01

51079-0299-63

51079-0299-66

51079-0299-17

51079-0299-20


Excipients: Ammonium Chloride ; Silicon Dioxide, Colloidal ; Croscarmellose Sodium (Type A) ; Magnesium Stearate ; Cellulose, Microcrystalline ; Sodium Lauryl Sulfate

Regulatory Status: RX

Availability: United States

Product ID 5150133

Contact: Mylan Pharmaceuticals



Copyright © 1974-2011 Thomson Reuters. All rights reserved.

DRUG POWDER ANALYSIS FORM

SAMPLE # [REDACTED] AGENCY Boston ANALYST ASD

No. of samples tested: Evidence Wt.

PHYSICAL DESCRIPTION:

3 white oval tablets
in bag

Gross Wt ():

Gross Wt ():

Pkg. Wt:

Net Wt: 4.8404

84443

600

PRELIMINARY TESTS

Spot Tests

Cobalt
Thiocyanate ()

Marquis

Froehde's

Mecke's

Microcrystalline Tests

Gold
Chloride

TLTA ()

OTHER TESTS

Monomer A/P

PRELIMINARY TEST RESULTS

RESULTS gabapentin

DATE 08-31-11

GC/MS CONFIRMATORY TEST

RESULTS gabapentin

MS
OPERATOR DEM

DATE 9-19-11

Drug Information

DRUG NAME: GABAPENTIN

INGREDIENTS: GABAPENTIN -- 600 MG

Color: White

Shape: Oval

Imprint: Hourglass 4443 600
Imprinted with a hourglass logo and 4443 on one side and 600 on the reverse side. , 4443 Imprinted with a hourglass logo and 4443 on one side and 600 on the reverse side. , 600

Form: Oral Tablet

Available-Container-Size: Bottle of 100 , Strip of 100 , Bottle of 1000 , Bottle of 500

AAPCC Code: 077730 - Other Types of Anticonvulsant (Excluding Barbiturates)

NDC: 0093-4443-01

0093-4443-05

0172-4443-60

0172-4443-10

0172-4443-80

0172-4443-70

0093-4443-10

Excipients: Lactose (Monohydrate) ; Cellulose, Microcrystalline ; Polyethylene Glycol

Regulatory Status: RX

Availability: United States

Product ID 6738342

Contact: Teva Pharmaceuticals

Copyright © 1974-2011 Thomson Reuters. All rights reserved.

DRUG POWDER ANALYSIS FORM

SAMPLE # AGENCY Boston ANALYST AS?

No. of samples tested: Evidence Wt.

PHYSICAL DESCRIPTION:


4 white round tablet
17 ipb

Gross Wt ():

Gross Wt ():

Pkg. Wt:

Net Wt: 0.3679



PRELIMINARY TESTS

Spot Tests

Cobalt
Thiocyanate ()

Marquis

Froehde's

Mecke's

Microcrystalline Tests

Gold
Chloride

TLTA ()

OTHER TESTS

Muonidesc A/P

PRELIMINARY TEST RESULTS

RESULTS clonidine

DATE 08-31-11

GC/MS CONFIRMATORY TEST

RESULTS clonidine

MS
OPERATOR nem

DATE 9-19-11

Drug Information

DRUG NAME: CLONIDINE HYDROCHLORIDE

INGREDIENTS: CLONIDINE HYDROCHLORIDE -- 0.1 MG

Color: White

Shape: Circle

Imprint: MYLAN 152

Form: Oral Tablet

Available-Container-Size: Strip of 90 , Strip of 25 , Bottle of 100 , Bottle of 1000 , Strip of 1 , Strip of 30 , Strip of 100

AAPCC Code: 201081 - Clonidine

NDC: 51079-0299-62

51079-0299-19

00378-0152-01

00378-0152-10

51079-0299-01

51079-0299-63

51079-0299-66

51079-0299-17

51079-0299-20


Excipients: Ammonium Chloride ; Silicon Dioxide, Colloidal ; Croscarmellose Sodium (Type A) ; Magnesium Stearate ; Cellulose, Microcrystalline ; Sodium Lauryl Sulfate

Regulatory Status: RX

Availability: United States

Product ID 5150133

Contact: Mylan Pharmaceuticals



Copyright © 1974-2011 Thomson Reuters. All rights reserved.

DRUG POWDER ANALYSIS FORM

SAMPLE # [REDACTED] AGENCY Rosier ANALYST AS9

No. of samples tested: Evidence Wt.

PHYSICAL DESCRIPTION:

39 white oval tablet
in pb

600

84443

Gross Wt ():

Gross Wt ():

Pkg. Wt:

Net Wt: 37.9467

PRELIMINARY TESTS

Spot Tests

Cobalt
Thiocyanate ()

Marquis

Froehde's

Mecke's

Microcrystalline Tests

Gold
Chloride

TLTA ()

OTHER TESTS

Monomelic AIP

PRELIMINARY TEST RESULTS

RESULTS gabapentin

DATE 08-31-11

GC/MS CONFIRMATORY TEST

RESULTS gabapentin

MS
OPERATOR Nem

DATE 9-19-11

Drug Information

DRUG NAME: GABAPENTIN

INGREDIENTS: GABAPENTIN -- 600 MG

Color: White

Shape: Oval

Imprint: Hourglass 4443 600
Imprinted with a hourglass logo and 4443 on one side and 600 on the reverse side. , 4443 Imprinted with a hourglass logo and 4443 on one side and 600 on the reverse side. , 600

Form: Oral Tablet

Available-Container-Size: Bottle of 100 , Strip of 100 , Bottle of 1000 , Bottle of 500

AAPCC Code: 077730 - Other Types of Anticonvulsant (Excluding Barbiturates)

NDC: 0093-4443-01

0093-4443-05

0172-4443-60

0172-4443-10

0172-4443-80

0172-4443-70

0093-4443-10

Excipients: Lactose (Monohydrate) ; Cellulose, Microcrystalline ; Polyethylene Glycol

Regulatory Status: RX

Availability: United States

Product ID 6738342

Contact: Teva Pharmaceuticals

Copyright © 1974-2011 Thomson Reuters. All rights reserved.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName DataFile LimsID	Method	Inj	SampleType	InjVolume
----	-----	-----	-----	-----	-----	-----
1	Vial 1	Blank	GENSCAN	1	Sample	1.0
2	Vial 2	Coke/Codeine Std	GENSCAN	1	Sample	1.0
3	Vial 3	Blank	GENSCAN	1	Sample	1.0
4	Vial 4	██████████	GENSCAN	1	Sample	1.0
5	Vial 5	Blank	GENSCAN	1	Sample	1.0
6	Vial 6	██████████	GENSCAN	1	Sample	1.0
7	Vial 7	Blank	GENSCAN	1	Sample	1.0
8	Vial 8	██████████	GENSCAN	1	Sample	1.0
9	Vial 9	Blank	GENSCAN	1	Sample	1.0
10	Vial 10	██████████	GENSCAN	1	Sample	1.0
11	Vial 11	Blank	GENSCAN	1	Sample	1.0


Sequence Table (Back Injector):

No entries - empty table!

Sample Name: Blank

```
=====
Acq. Operator   : ASD                      Seq. Line :    1
Acq. Instrument : DrugLab GC #2            Location  : Vial 1
Injection Date  : 8/31/2011 11:14:44 AM    Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.72019e4	5.22043e4	1.000e2

Totals : 5.72019e4 5.22043e4

```
=====
*** End of Report ***
```



Sample Name: Coke/Codeine Std

```

=====
Acq. Operator   : ASD                      Seq. Line :    2
Acq. Instrument : DrugLab GC #2            Location  : Vial 2
Injection Date  : 8/31/2011 11:49:27 AM    Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====

```



```

=====
                          Area Percent Report
=====

```

```

Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:      :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.77112e4	5.33250e4	98.66110
2	6.921	1	BB	403.30124	531.31085	0.68947
3	7.471	1	BB	379.88098	479.69025	0.64943

```
Totals :                5.84944e4  5.43360e4
```

```


=====
*** End of Report ***

```

Sample Name: Blank

```
=====
Acq. Operator   : ASD                      Seq. Line :    3
Acq. Instrument : DrugLab GC #2            Location  : Vial 3
Injection Date  : 8/31/2011 12:24:17 PM    Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:      :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.89751e4	5.57284e4	1.000e2

Totals : 5.89751e4 5.57284e4

```
=====
                        *** End of Report ***
=====
```


Sample Name: XXXXXXXXXX

```

=====
Acq. Operator   : ASD                      Seq. Line :    4
Acq. Instrument : DrugLab GC #2            Location  : Vial 4
Injection Date  : 8/31/2011 12:59:06 PM    Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====

```



```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.065	1	BB S	4.62303e4	4.87146e4	99.71723
2	6.000	1	BB	54.34018	47.75171	0.11721
3	6.612	1	VB	76.75762	74.00050	0.16556

```
Totals :                4.63614e4  4.88363e4
```

```


=====
*** End of Report ***

```

Sample Name: Blank

```
=====
Acq. Operator   : ASD                      Seq. Line :    5
Acq. Instrument : DrugLab GC #2             Location  : Vial 5
Injection Date  : 8/31/2011 1:33:55 PM      Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.98988e4	5.68185e4	1.000e2


Totals : 5.98988e4 5.68185e4

```
=====
*** End of Report ***
=====
```

Sample Name: XXXXXXXXXX

```
=====
Acq. Operator   : ASD                      Seq. Line :    6
Acq. Instrument : DrugLab GC #2            Location  : Vial 6
Injection Date  : 8/31/2011 2:08:45 PM      Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                          Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.63668e4	5.57624e4	99.76234
2	6.000	1	BB	51.89137	41.00683	0.09184
3	6.612	1	VB	82.39043	83.23665	0.14582


Totals : 5.65011e4 5.58867e4

```
=====
*** End of Report ***
```

Sample Name: Blank

```
=====
Acq. Operator   : ASD                      Seq. Line :    7
Acq. Instrument : DrugLab GC #2            Location  : Vial 7
Injection Date  : 8/31/2011 2:43:35 PM      Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.063	1	BB S	6.31883e4	5.75707e4	1.000e2

Totals : 6.31883e4 5.75707e4

```
=====
*** End of Report ***
```


Sample Name: XXXXXXXXXX

```

=====
Acq. Operator   : ASD                      Seq. Line :    8
Acq. Instrument : DrugLab GC #2            Location  : Vial 8
Injection Date  : 8/31/2011 3:18:24 PM      Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====

```



```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:      :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.063	1	BB S	5.81392e4	5.89850e4	94.68585
2	4.797	1	BB	2998.32910	1039.57239	4.88309
3	4.931	1	BB	126.71463	28.22244	0.20637
4	7.947	1	VB	74.62600	83.73415	0.12154
5	8.010	1	BB	28.25889	25.39233	0.04602
6	10.499	1	BB	35.08017	16.88823	0.05713

```
Totals :                6.14022e4  6.01788e4
```

```


=====
*** End of Report ***
=====

```

Sample Name: Blank

```
=====
Acq. Operator   : ASD                      Seq. Line :    9
Acq. Instrument : DrugLab GC #2             Location  : Vial 9
Injection Date  : 8/31/2011 3:53:13 PM      Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                        Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.062	1	BB S	5.46462e4	5.82483e4	1.000e2


Totals : 5.46462e4 5.82483e4

```
=====
*** End of Report ***
=====
```


Sample Name: XXXXXXXXXX

```
=====
Acq. Operator   : ASD                      Seq. Line :   10
Acq. Instrument : DrugLab GC #2            Location  : Vial 10
Injection Date  : 8/31/2011 4:28:03 PM     Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====
```



```
=====
                          Area Percent Report
=====
```

```
Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.063	1	BB S	5.86294e4	5.89250e4	96.91959
2	4.777	1	BB	1768.78735	786.53101	2.92396
3	7.947	1	VB	55.12289	58.79731	0.09112
4	8.011	1	BB	16.27322	15.07630	0.02690
5	10.500	1	BB	23.24638	10.56067	0.03843

Totals : 6.04928e4 5.97959e4

```
=====
*** End of Report ***
```


Sample Name: Blank

```

=====
Acq. Operator   : ASD                      Seq. Line :   11
Acq. Instrument : DrugLab GC #2           Location  : Vial 11
Injection Date  : 8/31/2011 5:02:53 PM    Inj       :    1
                                           Inj Volume: 1 µl

Sequence File   : C:\CHEM32\1\SEQUENCE\DEFAULT.S
Method          : C:\CHEM32\1\METHODS\GENSCAN.M
Last changed    : 3/9/2011 1:49:11 PM by ASD
=====

```



```

=====
                        Area Percent Report
=====

```

```

Sorted By           :      Retention Time
Multiplier:         :      1.0000
Dilution:           :      1.0000
Sample Amount:       :      1.00000 [ng/ul]   (not used in calc.)
Do not use Multiplier & Dilution Factor with ISTDs

```

Signal 1: FID1 A,

Peak #	RetTime [min]	Sig	Type	Area [pA*s]	Height [pA]	Area %
1	1.063	1	BB S	5.77777e4	5.95097e4	1.000e2

Totals : 5.77777e4 5.95097e4

```


=====
                        *** End of Report ***
=====

```

✓
NEM
01-20-12

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801074.D
Operator : NEM
Date Acquired : 15 Sep 2011 6:58
Sample Name : BLANK
Submitted by :
Vial Number : 2
AcquisitionMeth: GENSCAN.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:


File Name : F:\Q3-2011\System7\09_14_11\801075.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 7:22
 Sample Name : CLONIDINE STD PREP 9-9-2011
 Submitted by :
 Vial Number : 75
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE



Ret. Time	Area	Area %	Ratio %
7.363	2592278	100.00	100.00

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801077.D
Operator : NEM
Date Acquired : 15 Sep 2011 8:10
Sample Name : BLANK
Submitted by : ASD
Vial Number : 2
AcquisitionMeth: GENSCAN.M
Integrator : RTE



Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801075.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 7:22
 Sample Name : CLONIDINE STD PREP 9-9-2011
 Submitted by :
 Vial Number : 75
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE

Search Libraries: C:\Database\SLI.L Minimum Quality: 80
 C:\Database\NIST05a.L Minimum Quality: 80
 C:\Database\PMW_TOX2.L


PK#	RT	Library/ID	CAS#	Qual
1	7.36	C:\Database\NIST05a.L		
		Clonidine	004205-90-7	98
		Clonidine	004205-90-7	94
		Clonidine	004205-90-7	91



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801078.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 8:33
 Sample Name : XXXXXXXXXX
 Submitted by : ASD
 Vial Number : 78
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE



Ret. Time	Area	Area %	Ratio %
3.612	112535	4.49	4.70
7.369	2392946	95.51	100.00



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801078.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 8:33
 Sample Name : XXXXXXXXXX
 Submitted by : ASD
 Vial Number : 78
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE

Search Libraries: C:\Database\SLI.L Minimum Quality: 80
 C:\Database\NIST05a.L Minimum Quality: 80
 C:\Database\PMW_TOX2.L

PK#	RT	Library/ID	CAS#	Qual
1	3.61	C:\Database\NIST05a.L		
		Benzenecarboxylic acid	000065-85-0	91
		Benzenecarboxylic acid	000065-85-0	91
		Heptanediamide, N,N'-di-benzoyloxy-	1000253-26-4	83





Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801078.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 8:33
 Sample Name : XXXXXXXXXX
 Submitted by : ASD
 Vial Number : 78
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE

Search Libraries: C:\Database\SLI.L Minimum Quality: 80
 C:\Database\NIST05a.L Minimum Quality: 80
 C:\Database\PMW_TOX2.L


PK#	RT	Library/ID	CAS#	Qual
2	7.37	C:\Database\NIST05a.L		
		Clonidine	004205-90-7	98
		Clonidine	004205-90-7	95
		Clonidine	004205-90-7	78



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801080.D
Operator : NEM
Date Acquired : 15 Sep 2011 9:21
Sample Name : BLANK
Submitted by : ASD
Vial Number : 2
AcquisitionMeth: GENSCAN.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801081.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 9:45
 Sample Name : XXXXXXXXXX
 Submitted by : ASD
 Vial Number : 81
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE



Ret. Time	Area	Area %	Ratio %
7.363	2752900	100.00	100.00



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801081.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 9:45
 Sample Name : XXXXXXXXXX
 Submitted by : ASD
 Vial Number : 81
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE

Search Libraries: C:\Database\SLI.L Minimum Quality: 80
 C:\Database\NIST05a.L Minimum Quality: 80
 C:\Database\PMW_TOX2.L


PK#	RT	Library/ID	CAS#	Qual
1	7.36	C:\Database\NIST05a.L		
		Clonidine	004205-90-7	98
		Clonidine	004205-90-7	91
		Clonidine	004205-90-7	78



Area Percent / Library Search Report

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801083.D
 Operator : NEM
 Date Acquired : 15 Sep 2011 10:33
 Sample Name : BLANK
 Submitted by :
 Vial Number : 2
 AcquisitionMeth: GENSCAN.M
 Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801084.D
Operator : NEM
Date Acquired : 15 Sep 2011 10:57
Sample Name : CLONIDINE STD PREP 9-9-2011
Submitted by :
Vial Number : 75
AcquisitionMeth: GENSCAN.M
Integrator : RTE





Ret. Time	Area	Area %	Ratio %
7.363	2312206	100.00	100.00

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801084.D
Operator : NEM
Date Acquired : 15 Sep 2011 10:57
Sample Name : CLONIDINE STD PREP 9-9-2011
Submitted by :
Vial Number : 75
AcquisitionMeth: GENSCAN.M
Integrator : RTE


Search Libraries: C:\Database\SLI.L Minimum Quality: 80
C:\Database\NIST05a.L Minimum Quality: 80
C:\Database\PMW_TOX2.L

PK#	RT	Library/ID	CAS#	Qual
1	7.36	C:\Database\NIST05a.L		
		Clonidine	004205-90-7	98
		Clonidine	004205-90-7	90
		Clonidine	004205-90-7	78



Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801086.D
Operator : NEM
Date Acquired : 15 Sep 2011 11:43
Sample Name : BLANK
Submitted by :
Vial Number : 2
AcquisitionMeth: MDMA.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801087.D
Operator : NEM
Date Acquired : 15 Sep 2011 13:11
Sample Name : GABAPENTIN STD
Submitted by :
Vial Number : 87
AcquisitionMeth: MDMA.M
Integrator : RTE





Ret. Time	Area	Area %	Ratio %
2.465	1106473	100.00	100.00

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801087.D
Operator : NEM
Date Acquired : 15 Sep 2011 13:11
Sample Name : GABAPENTIN STD
Submitted by :
Vial Number : 87
AcquisitionMeth: MDMA.M
Integrator : RTE


Search Libraries: C:\Database\SLI.L Minimum Quality: 80
C:\Database\NIST05a.L Minimum Quality: 80
C:\Database\PMW_TOX2.L

PK#	RT	Library/ID	CAS#	Qual
1	2.47	C:\Database\PMW_TOX2.L		
		Gabapentin -H2O	000000-00-0	97
		Lenacil	002164-08-1	12
		Chlormezanone	000080-77-3	2



Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801089.D
Operator : NEM
Date Acquired : 15 Sep 2011 13:34
Sample Name : BLANK
Submitted by : ASD
Vial Number : 2
AcquisitionMeth: MDMA.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801090.D
Operator : NEM
Date Acquired : 15 Sep 2011 13:45
Sample Name : XXXXXXXXXX
Submitted by : ASD
Vial Number : 90
AcquisitionMeth: MDMA.M
Integrator : RTE





Ret. Time	Area	Area %	Ratio %
2.471	4191550	100.00	100.00

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801090.D
Operator : NEM
Date Acquired : 15 Sep 2011 13:45
Sample Name : XXXXXXXXXX
Submitted by : ASD
Vial Number : 90
AcquisitionMeth: MDMA.M
Integrator : RTE


Search Libraries: C:\Database\SLI.L Minimum Quality: 80
C:\Database\NIST05a.L Minimum Quality: 80
C:\Database\PMW_TOX2.L

PK#	RT	Library/ID	CAS#	Qual
1	2.47	C:\Database\PMW_TOX2.L		
		Gabapentin -H2O	000000-00-0	97
		Lenacil	002164-08-1	12
		Chlormezanone artifact	000000-00-0	1



Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801092.D
Operator : NEM
Date Acquired : 15 Sep 2011 14:08
Sample Name : BLANK
Submitted by : ASD
Vial Number : 2
AcquisitionMeth: MDMA.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801095.D
Operator : NEM
Date Acquired : 15 Sep 2011 14:43
Sample Name : BLANK
Submitted by :
Vial Number : 2
AcquisitionMeth: MDMA.M
Integrator : RTE




Ret. Time	Area	Area %	Ratio %
-----------	------	--------	---------

NO INTEGRATED PEAKS

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801096.D
Operator : NEM
Date Acquired : 15 Sep 2011 14:54
Sample Name : GABAPENTIN STD
Submitted by :
Vial Number : 87
AcquisitionMeth: MDMA.M
Integrator : RTE



Ret. Time	Area	Area %	Ratio %
2.460	1201929	100.00	100.00

Information from Data File:

File Name : F:\Q3-2011\System7\09_14_11\801096.D
Operator : NEM
Date Acquired : 15 Sep 2011 14:54
Sample Name : GABAPENTIN STD
Submitted by :
Vial Number : 87
AcquisitionMeth: MDMA.M
Integrator : RTE

Search Libraries: C:\Database\SLI.L Minimum Quality: 80
C:\Database\NIST05a.L Minimum Quality: 80
C:\Database\PMW_TOX2.L

PK#	RT	Library/ID	CAS#	Qual
1	2.46	C:\Database\PMW_TOX2.L		
		Gabapentin -H2O	000000-00-0	96
		Lenacil	002164-08-1	8

